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Instruction For Use

valid for all products **Urological Catheters (non-sterile/BULK)** and **Catheter Plugs (non-sterile/BULK)**

1. Manufacturer: MEDINORM Medizintechnik GmbH, Gewerbepark 7-9, 66583 Spiesen-Elversberg, Germany
2. **Catheter**; for urologic use; Accessory: Catheter Plugs
3. Intended use of sterile device: Urine removal from bladder. Catheter stopper (accessory) for temporary hygienic closure. Indication: Necessary bladder removal, Urine recovery or accounting. Contraindication: Urinary inflammation, discontinuity or strictur of Urethra. Specification of patient target group not applicable.
4. **Non-Sterile Medical Device**; Unpack under Cleanroom condition; Do not use by damaged inner packaging! For single use in sterile version!
5. Residual risk or risk by unintended use: risk of infection!
6. Protect from strong light!
7. Unused product is recyclable; used product is to be disposed contamination-protected!
8. In case of incidence consider national notification requirements!
9. Use through professional persons only

Note: Intended use only possible STERILE!

Note: Necessity of existing validated sterile packaging process!

Note: Sterilisation (EO) through validated process necessary!

Note: Sterilisation once possible under use of below instructed parameters!

Note: Repeated or Re-Sterilisation must be validated separately!

Process parameter to be used:

Preconditioning:	
Chamber temperature at end of conditioning phase(°C)	40 ±5
Relative humidity in chamber end of conditioning phase (rH%)	60 ±15
Preconditioning time (h)	6 – 48
Sterilisation:	
Initial vacuum (mbar _{abs})	100 +0/-30
Relative humidity in chamber (at the end of conditioning) (rH%)	70 ±20
Pressure increase after EO injection (mbar _{abs})	390 +30/-0
Gas exposure time (min)	240 +5/-0
Chamber temperature during gas exposure time (°C)	40 +6/-0
Post conditioning:	
Temperature post conditioning chamber (°C)	39 ± 6
Time of desorption (h)	≥48

Explanation of symbols used: <https://medinorm.eu/IFU/Symbols.pdf>

This user manual can be requested in paper form from the manufacturer MEDINORM free of charge. Delivery takes place within seven calendar days. Always have a current version in an appropriate form in your emergency management system!