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Instruction For Use

valid for all products **High and Low Vacuum Wound Drainage Systems (non-sterile/BULK)**

1. Manufacturer: MEDINORM Medizintechnik GmbH, Gewerbepark 7-9, 66583 Spiesen-Elversberg, Germany
2. **High / Low Vacuum Wound Drainage System**; with Drain-connector; for acute wound healing;
3. Intended use of sterile device: Secretion derive from wound field by vacuum. Indication: Initiative wound healing. Contraindication: disordered blood clotting, Deformity, Organ contact. Specification of patient target group only regarding to Medical Result.
4. **Non-Sterile Medical Device**; Unpack under Cleanroom condition; Do not use by damaged inner packaging! For single use in sterile version!
5. Residual risk or risk by unintended use: risk of infection!
6. Protect from strong light!
7. Unused product is recyclable; used product is to be disposed contamination-protected!
8. Risk of infection! Risk of stitch damage! Observe National Pinprick Guidelines and Regulations!
9. In case of incidence consider national notification requirements!
10. Use through professional persons only:
 - I. Choose intended drainage tube.
 - II. Place the drainage, take care of airtight wound closure.
 - III. Unpack Drainage bottle (bellow contracted!) and connect tube to bottle; if already prepared assure tight fitting.
 - IV. Cut the Drainage connector corresponding to chosen Drainage size.
 - V. Insert the Drainage into Connector, take care of tight fitting.
 - VI. Open distal clamp (bottle side) = System activated!
 - VII. Before any necessary manipulation close proximal clamp (tube side), then distal clamp (bottle side).
 - VIII. Appropriate waste disposal according to national requirements.Bottle exchange: See VII. + change bottle + See VI + open proximal clamp.

Note: Intended use only possible STERILE!

Note: Necessity of existing validated sterile packaging process!

Note: Sterilisation (EO) through validated process necessary!

Note: Sterilisation once possible under use of below instructed parameters!

Note: Repeated or Re-Sterilisation must be validated separately!

Process parameter to be used:

Preconditioning:	
Chamber temperature at end of conditioning phase(°C)	40 ±5
Relative humidity in chamber end of conditioning phase (rH%)	60 ±15
Preconditioning time (h)	6 – 48
Sterilisation:	
Initial vacuum (mbar _{abs})	100 +0/-30
Relative humidity in chamber (at the end of conditioning) (rH%)	70 ± 20
Pressure increase after EO injection (mbar _{abs})	460 ± 10
Gas exposure time (min)	360 + 5/-0
Chamber temperature during gas exposure time (°C)	40 +5/-0
Post conditioning:	
Temperature post conditioning chamber (°C)	40 ± 5
Time of desorption (h)	≥48

Explanation of symbols used: <https://medinorm.eu/IFU/Symbols.pdf>

This user manual can be requested in paper form from the manufacturer MEDINORM free of charge. Delivery takes place within seven calendar days. Always have a current version in an appropriate form in your emergency management system!